



# CROWNTEC | Compatibility Overview

3D PRINTER	CLEANING	POST-CURING*	POST-TREATMENT
 ACKURETTA SOL 3 DENTIQ 3	<p><b>CLEANING BY HAND</b></p> <p>Clean the print object with brushes and cloths minimal use of IPA (96%) from excess resin.</p> <p>Air dry the surface and inside of the print object.</p>	<b>1</b> OtoFlash (NK-Optik) 2 x 2000 Flashes	<p><b>RECOMMENDED WATER BATH**</b></p> <p>Place the print object in boiling water (100°C / 212°F) for 2 minutes after post-curing to finalize the color.</p> 
 ASIGA MAX UV 1 2 PRO 4K 1 2		<b>2</b> HiLite Power (Kulzer) 2 x 180 s	
 DentaFab		<b>3</b> Curie (Ackuretta) 2 x 3 min   T:6 P:10 D:16 B:ON	
 Desktop Health		<b>4</b> Cure (DentaFab) 1 x 10 min	
 DMG 3Delite 1 2 9 3Demax 1 2 9		<b>5</b> Photopol (Dentalfarm) 2 x 4 min with vacuum	
 MIICRAFT Prime 1		<b>6</b> BB Compact (Meccatronicore) 1 x 6 min   Power 100%   Temp. 20°C / 68°F Plate on top rail	<p><b>OPTIONAL POLIMERIZATION LAMP**</b></p> <p>Expose each side of the print object for 2 x 20 seconds at full power.</p> 
 Next Dent by 3D SYSTEMS NextDent 5100 7		<b>7</b> LC-3DPrint Box (NextDent) 1 x 30 min	
 phrozen Sonic 4K 2022 8 Sonic 4K XL 2022 8		<b>8</b> Cure (Phrozen) 2 x 5 min	
 rapidshape D20 II, D30 II, D40 II, D10+, D20+ cartridge, D20+, D30+, D50+ 1 2 9		<b>9</b> RS cure (Rapid Shape) 1 x 6 min	
 SHINING 3D		<b>10</b> Fab Cure (Shining) 2 x 10 min	

\*Recommended polymerization devices for light-curing materials such as OtoFlash and HiLite Power reach a wavelength range of 320 - 500nm. Other polymerization devices mentioned above do not reach the upper wavelength range and do not completely finish the esthetic color finalization process (no influence on the physical properties of the material).

\*\*To accelerate the color finalization, it is recommended to place the print object in boiling water (100°C / 212°F) for 2 minutes. A polymerization lamp such as a Bluephase® G2 from Ivoclar Vivadent can be used optionally (2 x 20 s full power on the print object per side) to cover the higher wavelength ranges (up to 500nm).

Please note: This compatibility overview does not replace the instructions for use. Please read the instructions for use carefully.

Edited: 09/2022 | D600243 | Class IIa Medical Device

More information at [saremco.ch/en/downloads/instruction-for-use/](https://saremco.ch/en/downloads/instruction-for-use/)